



## Total RNA extraction for turbot eggs (DevMap)

This protocol is a modification of the original NMBU protocol "Total RNA extraction for tissues (BodyMap)"

[https://data.faang.org/api/fire\\_api/experiments/NMBU\\_SOP\\_RNAextraction\\_protocol\\_20200503.pdf](https://data.faang.org/api/fire_api/experiments/NMBU_SOP_RNAextraction_protocol_20200503.pdf)

### Preparation

- Pre-chill a fixed rotor centrifuge with capacity for 2 mL tubes to 4°C, under the fume hood
- All steps should be done at room temperature, except the 1<sup>st</sup> centrifugation (step 6)
- All steps should be done under the fume hood
- Prior to first use, prepare the stock DNase I mix, the buffer RWT and RPE according to manufacturer's instructions

### Protocol

1. Remove all TRIzol from the 2 mL safe-lock tube containing the eggs. Add 700 ul of Qiazol Lysis Reagent and 2 zirconia beads 2.8-3.3 mm.
2. Under the fume hood, place the tube in TissueLyser II and disrupt for 2 minutes at 20 Hz. Let the samples rest for 1 minute and disrupt again for 2 minutes at 20 Hz.
3. Let the tubes stand vertically for 5 minutes.
4. Add 140 ul of chloroform and vortex for 15 s.
5. Let the tubes stand vertically for 2-3 minutes.
6. Centrifuge 15 minutes 12,000 x g at 4°C. After this step, heat the centrifuge up to room temperature.
7. Carefully remove the tube from centrifuge and transfer the upper aqueous phase (approximately 350 ul) to a new 1.5 mL tube.
8. Add exactly 1.5x volume of 100% ethanol (usually 525 ul), mix well by pipetting and continue with the next step without delay.
9. Pipet up to 700 ul of the sample into a RNeasy Mini spin column placed into a 2 mL collection tube. Centrifuge at >8000 x g for 15 s at room temperature and discard the flow-through.
10. Repeat the previous step with the remaining volume of sample and with the same column. Discard the flow-through.
11. Add 350 ul of buffer RWT in the column and centrifuge for 15 s at > 8000 x g. Discard the flow-through.
12. Add 80 ul of DNase I mix and incubate it at room temperature for 15 min.
13. Pipet 350 ul of buffer RWT in the column and centrifuge for 15 s at > 8000 x g. Discard the flow-through.
14. Pipet 500 ul of buffer RPE in the column and centrifuge for 15 s at > 8000 x g. Discard the flow-through.

15. Pipet another 500 ul of buffer RPE in the column and centrifuge for **2 minutes** at > 8000 x g. Discard the flow-through.
16. Place the column in a new 2 ml collection tube and centrifuge 1 minute at > 13000 x g, to remove any remains of previous buffers.
17. Place the column in a new 1.5 ml tube. Add **50** ul of RNase free water and centrifuge 1 minute at > 8000 g. Discard the column and place the tube with the flow-through on ice.
18. Quantify your RNA on Nanodrop and check its RIN value and profile on Bioanalyzer. Freeze the remaining 50 ul at - 80 °C as soon as you can.

## Reagents and Equipment

Reagent	Reference
<b>miRNeasy Mini kit (Qiagen)</b>	217004
<b>RNase free DNase set (Qiagen)</b>	79254
<b>2.8-3.3mm zirconia beads</b>	9738463
<b>Bioanalyzer RNA 6000 Nano kit</b>	5067-1511
TissueLyser II	
Microcentrifuge thermoregulated	
<b>Ethanol 100%</b>	
<b>Chloroform</b>	
<b>Ice</b>	